

REBEL

COGENT AMERICA INC

Suite 104, 31007 IH-10, West Boerne,
Kendall County, Texas, 78006

T. 830-368-4210

E. hello@cogentamerica.com

www.cogentamerica.com



WHAT MAKES REBEL SO SPECIAL?



COMPACT, INTEGRATED DESIGN

- Low part count, low complexity, low weight, low cost
- Specifically designed for woman & teen
- Hand width 71.4mm, 160mm length, knuckle diameter 19mm

MARKET LEADING AESTHETICS

- Designed to be a highly desirable product, not a medical device

WATERPROOF

- No glove required, no user worries, practical and easy to clean

HIGHLY FUNCTIONAL

- Passive flex wrist to reduce shoulder compensatory movement.
- Tap to move thumb position, quick function for all users.
- Passive flex digits to protect from impact damage.
- Force vectoring grip optimisation, low cognitive load gripping small items

